

Introduction to RNA-Seq – Quality Control

Wandrille Duchemin


















"Raw data": FASTQ format

Paired data: two separate files for forward and reverse with same ordering


"Raw data": FASTQ format

Paired data: two separate files for forward and reverse with same ordering

"Raw data": FASTQ format - header



"Raw data": FASTQ format - header



Depends on sequencing technology used, changed several times by Illumina and others

TCTCNAGATAAAATCAAAACCAACAGAGAGTCTAGAATAAAAGTGAATAG

Each nucleotide is associated to a quality value

"Raw data": FASTQ format - Phred score

- ASCII encoded
- Represents probability (p) that base call is incorrect
 - Quality (Q) = -10 $\log_{10}(p)$

P-value	Phred	Probability of incorrect base call	Base call accuracy		
1E-01	10	1 in 10	90%		
1E-02	20	1 in 100	99%		
1E-03	30	1 in 1000	99.9%		
1E-04	40	1 in 10,000	99.99%		

TCTCNAGATAAAATCAAAACCAACAGAGAGTCTAGAATAAAAGTGAATAG

@BF#2ADHHHHHJJJJIJJJJJJJJJJJJJHIIGIHIIIIJJHIHIJJJ

Q	ASCI	I P	Q	ASCII	Р	Q	ASCII	Р	Q	ASCII	P
1	"	0.79433	12	-	0.06310	23	8	0.00501	34	С	0.00040
2	#	0.63096	13		0.05012	24	9	0.00398	35	D	0.00032
3	\$	0.50119	14	/	0.03981	25	:	0.00316	36	E	0.00025
4	%	0.39811	15	0	0.03162	26	;	0.00251	37	F	0.00020
5	&	0.31623	16	1	0.02512	27	<	0.00200	38	G	0.00016
6	•	0.25119	17	2	0.01995	28	=	0.00158	39	H	0.00013
7	(0.19953	18	3	0.01585	29	>	0.00126	40	I	0.00010
8)	0.15849	19	4	0.01259	30		0.00100	41	J	0.00008
9	*	0.12589	20	5	0.01000	31	@	0.00079			
10	+	0.10000	21	6	0.00794	32	А	0.00063			
10	l,	0.07943	22	7	0.00631	33	В	0.00050			

TCTCNAGATAAAATCAAAACCAACAGAGAGTCTAGAATAAAAGTGAATAG

Q	ASCII	Р	Ų	ASCII	Ρ	Ų	ASCII	Ρ	Q	ASCII	Ρ
1	"	0 79433	12	-	0.06310	23	8	0.00501	34	C	0.00040
2	#	0.63096	13		0.05012	24	9	0.00398	35	D	0.00032
- 5	>	0.50II9	14	/	0.03981	25	:	0.00316	36	E	0.00025
4	%	0.39811	15	0	0.03162	26	;	0.00251	37	F	0.00020
5	&	0.31623	16	1	0.02512	27	<	0.00200	38	G	0.00016
6		0.25119	17	2	0.01995	28	=	0.00158	39	H	0.00013
7	(0.19953	18	3	0.01585	29	>	0.00126	40	I	0.00010
8)	0.15849	19	4	0.01259	30	?	0.00100	41	J	0.00008
9	*	0.12589	20	5	0.01000	31	@	0.00079			
10	+	0.10000	21	6	0.00794	32	A	0.00063			
11	,	0.07943	22	7	0.00631	33	В	0.00050			

TCTCNAGATAAAATCAAAACCAACAGAGAGTCTAGAATAAAAGTGAATAG

Illumina v1.8 and later (ASCII_BASE=33)

Q	ASC	I P	Q	ASCII	P	Q	ASCII	Р	Q	ASCII	Р
1	"	0.79433	12	-	0.06310	23	8	0.00501	34	C	0.00040
2	#	0.63096	13		0.05012	24	9	0.00398	35	D	0.00032
3	\$	0.50119	14	/	0.03981	25	:	0.00316	36	E	0.00025
4	%	0.39811	15	0	0.03162	26	;	0.00251	37	F	0.00020
5	&	0.31623	16	1	0.02512	27	<	0.00200	38	G	0.00016
6		0.25119	17	2	0.01995	28	=	0.00158	30	Н	0.00013
7	(0.19953	18	3	0.01585	29	>	0.00126	40	I	0.00010
8)	0.15849	19	4	0.01259	30	;	0.00100	41	J	800000.0
9	*	0.12589	20	5	0.01000	31	@	0.00079			
10	9 +	0.10000	21	6	0.00794	32	A	0.00063			
1:	1,	0.07943	22	7	0.00631	33	В	0.00050			

"Raw data": FASTQ format – Phred +33/+64

Sanger, Illumina v1.3 to 1.7 (ASCII BASE=64)

Q	ASCII	P	Q	ASCII	Р	Q	ASCII	Р	Q	ASCII	Р
1	Α	0.79433	12	L	0.06310	23	W	0.00501	34	b	0.00040
2	В	0.63096	13	M	0.05012	24	X	0.00398	35	c	0.00032
3	C	0.50119	14	N	0.03981	25	Υ	0.00316	36	d	0.00025
4	D	0.39811	15	0	0.03162	26	Z	0.00251	37	e	0.00020
5	E	0.31623	16	P	0.02512	27	[0.00200	38	f	0.00016
6	F	0.25119	17	Q	0.01995	28	\	0.00158	39	g	0.00013
7	G	0.19953	18	R	0.01585	29]	0.00126	40	h	0.00010
8	H	0.15849	19	S	0.01259	30	^	0.00100			
9	I	0.12589	20	T	0.01000	31	_	0.00079			
10	J	0.10000	21	U	0.00794	32	~	0.00063			
11	K	0.07943	22	V	0.00631	33	a	0.00050			

Illumina v1.8 and later (ASCII_BASE=33)

Q	ASCII	Р									
1	"	0.79433	12	-	0.06310	23	8	0.00501	34	С	0.00040
2	#	0.63096	13		0.05012	24	9	0.00398	35	D	0.00032
3	\$	0.50119	14	/	0.03981	25	:	0.00316	36	E	0.00025
4	%	0.39811	15	0	0.03162	26	;	0.00251	37	F	0.00020
5	&	0.31623	16	1	0.02512	27	<	0.00200	38	G	0.00016
6	•	0.25119	17	2	0.01995	28	=	0.00158	39	H	0.00013
7	(0.19953	18	3	0.01585	29	>	0.00126	40	I	0.00010
8)	0.15849	19	4	0.01259	30	;	0.00100	41	J	0.00008
9	*	0.12589	20	5	0.01000	31	@	0.00079			
10	+	0.10000	21	6	0.00794	32	Α	0.00063			
11	,	0.07943	22	7	0.00631	33	В	0.00050			

Quality Control using FastQC

Help spot problems in the sequencer or in the starting library material

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/

- Generates an HTML report with various metrics
 - Average quality per position
 - GC% profile
 - Adapter presence
 - •
- input formats: FASTQ (gzip) , SAM, BAM

Quality Control using FastQC

Help spot problems in the sequencer or in the starting library material

https://www.bioinformatics.babraham.ac.uk/projects/fastqc/

- Generates an HTML report with various metrics
 - Average quality per position
 - GC% profile
 - Adapter presence
 - •
- input formats: FASTQ (gzip), SAM, BAM

Combining multiple reports : MultiQC

- FastQC: 1 report for each fastq file
- MultiQC combine individual reports into a single report

https://multiqc.info/


- MultiQC also works with other tools outputs:
 - Trimming output
 - Mapping output

•


PRACTICAL

Go to the website and do the QC practical

Per Base Sequence Quality




Per Tile Sequence Quality




Only present when the fastq id contains tile id


Per Sequence Quality Scores




Per Base Sequence Content




Per Sequence GC Content



Per Sequence GC Content





Duplicate Sequences



Overrepresented Sequences

Sequence	Count	Percentage	Possible Source
GATCGGAAGAGCACACGTCTGAACTCCAGTCACACAGTGATCTCGTATGC	355643	2.113348167370486	TruSeq Adapter, Index 5 (100% over 50bp)
${\tt AGATCGGAAGAGCACACGTCTGAACTCCAGTCACACAGTGATCTCGTATG}$	42318	0.2514675327414971	TruSeq Adapter, Index 5 (100% over 49bp)





Swiss Institute of Bioinformatics

Contributors:

Wandrille Duchemin Geoffrey Fucile Walid Gharib Pablo Escobar Lopez Mihaela Zavolan















